

SUPERNO₂VA® ET MASK

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE | MRI COMPATIBLE

The SuperNO₂VA® Et Mask is a non-invasive mask designed to deliver positive airway pressure and oxygen to patients at risk of upper airway obstruction during moderate to deep sedation procedures.



Adult Medium Mask
(SEM-20)

Adult Large Mask
(SEL-20)



PRE-OP
Supports
pre-oxygenation
to maintain O₂



INTRA-OP
Stents the airway open
during deep sedation
and allows for EtCO₂
monitoring



POST-OP
Provides support
during transport
and in PACU

- Designed to optimize oxygenation and ventilation
- Delivers nasal positive airway pressure and oxygen throughout the continuum of care
- Prevents unnecessary progression to invasive airway procedures
- Maintains open airway and minimizes upper airway obstruction
- Allows for full oral access, effectively share the airway

ITEM	DESCRIPTION	PK
SEM-20	SuperNO ₂ VA Et Mask, Adult Medium	20
SEL-20	SuperNO ₂ VA Et Mask, Adult Large	20
SNEM-20	SuperNO ₂ VA Et, Nasal PAP Ventilation Mask, Medium	20
SNEL-20	SuperNO ₂ VA Et, Nasal PAP Ventilation Mask, Large	20

SUPERNO₂VA® ET MASK AND SYSTEM

The SuperNO₂VA® Et Mask is also available packaged with a satellite system, including a hyperinflation system, supplemental O₂ tubing, and various gas sampling lines.

- Enables rescue ventilation without the need for capital equipment to ventilate the patient
- Offers gas sampling lines with male luer or reflective connectors to meet specific application needs
- Optimizes oxygenation and ventilation in both operating room and non-operating room anesthesia environments
- Prevents unnecessary progress to invasive airways

ITEM	DESCRIPTION	PK
SSEM-20	SuperNO ₂ VA Et System with Adult Medium Mask, 2 L Hyperinflation System, O ₂ Tubing and Male Gas Sampling Line	20
SSEL-20	SuperNO ₂ VA Et System with Adult Large Mask, 2 L Hyperinflation System, O ₂ Tubing, and Male Gas Sampling Line	20
SSNEM-20*	SuperNO ₂ VA Et System with Adult Medium Mask, 2 L Hyperinflation System, O ₂ Tubing, Male Gas Sampling Line and Reflective Gas Sampling Line	20
SSNEL-20*	SuperNO ₂ VA Et System with Adult Large Mask, 2 L Hyperinflation System, O ₂ tubing, Male Gas Sampling Line and Reflective Gas Sampling Line	20

*Available only in the United States



2L Hyperinflation Bag



Supplemental
O₂ Tubing



Reflective Gas
Sampling Line
(Compatible with
Microstream technology)



Gas Sampling
Line with Male
Connectors



Large or
Medium Mask