SunMed®


SunMed® is one of the world's largest and most comprehensive providers of anesthesia and respiratory products.

More than 30 trusted brand names reside under the SunMed umbrella.

SALTER LABS® Ethox medical ventlab INTUBRITE® Westmed Cookgas®

Global Manufacturing and Distribution Footprint

Low-cost manufacturing footprint in Asia and North America with China, US and EU distribution centers



6
Global Manufacturing Sites

820K+
Manufacturing Sq. Ft.

5
Global Distribution Sites

1
New Innovation Center

FDA Registered
All Manufacturing Sites

MDSAP Compliant
All Manufacturing Sites

ISO 13485
All Manufacturing Sites

Respiratory



NEBUTECH® NEBULIZER

- For ED, acute care, EMS and in-home therapy
- Optional mask, expiratory filter and PEP therapy
- Shorter treatment time helps promote patient compliance; 3 mL in less than 5 minutes at 8 LPM
- B/V filter option helps reduce the transmission of microbes and other particulate matter in the breathing system



HEART™ CONTINUOUS NEBULIZER

- One of the largest clinical coverage of standard and reflective connectors for $EtCO_2$ reducing health system costs and SKUs, while increasing patient care
- Comfort Soft Plus™ tubing drastically reduces HAPU while promoting patient comfort



VIXONE™ SMALL VOLUME JET NEBULIZER

- Popular acorn shape helps decrease medication waste by redirecting back to the jet
- Highly-efficient respirable mass output with MMAD of 2.7 microns promotes greater deposition
- Easy to fill cup with ¼ turn lock provides leak-proof seal
- Swivel mask feature allows treatment in upright positions or at angles up to 45 degrees
- 10 cc medication bowl capacity
- Very low dead volume (0.7 cc)

Masks



AEROSOL MASKS

- Swivel snout accommodates multiple patient positions
- Vented
- Regular and elongated models available
- Elongated models feature metal nose clips for secure fit

Pressure Infusion



INFU-SURG® AND INFU-SURG® CLARITY

- Safe, high quality, reliable device for intravenous and intra-arterial infusion treatment, including A-line pressure monitoring
- Pressure relief valve prevents overinflation (375 mm Hg pressure relief)
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- INFU-SURG CLARITY offers a transparent panel for clear visibility of solution bag

PS-10038 Rev2