SALTER LABS® PRO_2 ELITE PLUS™ OXYGEN ANALYZER

LATEX FREE

The Salter Labs® PRO_2 Elite Plus™ Oxygen Analyzer is a tool to be used by qualified personnel to spot-check oxygen concentration, flow, and pressure of a delivered air/oxygen mixture from an oxygen concentrator. The oxygen analyzer is not intended for use in continuous monitoring of oxygen delivery to patients.

FEATURES

- Ultrasound technology for oxygen concentration and flow measurements
- Silicon pressure sensor for pressure measurements
- Single gas inlet for all measurements
- No in-field calibration of oxygen/flow measurement required
- Auto-off after 5 minutes
- Self-diagnostics including error indication
- Battery powered with two (2) AA batteries (not included)
- Low battery indication
- Pressure measurements in PSI or kPa

INCLUDED COMPONENTS



- PRO_2 Elite Plus Oxygen Analyzer
- 2' Connection Tube with Trumpet and Universal Connectors
- Cord Clip



ITEM	DESCRIPTION	PK
PRO2PLUS-1	PRO_2 Elite Plus Oxygen Analyzer, 2' Connection Tube with Trumpet & Universal Connectors, & Cord Clip	1
PRO2PLUS-2	PRO_2 Elite Plus Oxygen Analyzer, 2' Connection Tube with Trumpet & Universal Connectors, & Cord Clip	1

PRO2PLUS-1 (U.S. version) *may be sold only in the United States.* ***PRO2PLUS-2 (International version)*** *is not intended for sale in the United States.*

LIT-20102-C1 Rev. 2